

































on-resonance will be higher than the fluorescence signal from an unpatterned glass substrate. Also, photobleaching is not an important issue for one-time scanning of DNA microarrays and protein microarrays.

## 6. Conclusions

This paper reports the design and fabrication of a PCEF surface that is fabricated upon a quartz substrate for low autofluorescence and high enhancement factors for simultaneous PC enhanced fluorophore excitation and PC enhanced fluorophore emission extraction. The PCEF surface gave a maximum enhancement factor of  $7500 \times$  for a  $\sim 50$  nm thick layer of LD-700 (concentration of 538 ng/ml) doped SU-8 layer. Using dose-response characterization of deposited PPL-Alexa 647 spots of variable tagged molecule concentration, a SNR improvement of  $330 \times$  on the PC was demonstrated for the concentration corresponding to the LOD on an unpatterned glass surface. The LOD on the PC slide was lowered by  $140 \times$  compared to the LOD of PPL-Alexa on the glass control.

The PCEF surface can be used to provide lower detection limits for broad classes of surface-based fluorescent assays for applications that include DNA microarrays for quantification of gene expression, protein microarrays for detection of disease biomarkers in blood, and next-generation DNA sequencing applications that utilize fluorescent tags.

## Acknowledgements

This work was supported by SRU Biosystems, the National Science Foundation (CBET 07-54122), and the National Institutes of Health (R01 GM086382). Any opinions, findings, and conclusions or recommendations expressed in this material are those of the authors and do not necessarily reflect the views of the National Science Foundation. The authors would like to thank Molecular Imprints, Inc. for their technical support. Contact \* B.T. Cunningham, Tel: + 1-217-265-6291; [bcunning@illinois.edu](mailto:bcunning@illinois.edu)